

## Additional equipment required for the Lab901 TapeStation®

It is imperative that only the consumables listed below are used in the TapeStation. Failure to do so may result in damage to the instrument.

The table below details the latest information regarding compatible tips, tubes and plates for the current configurations of the Lab901 TapeStation. All products listed below are available to purchase from Lab901 Ltd.

Equipment	Description	Catalogue Number	Quantity
Tips	ScreenTape® filtered loading tips	ST304	384 tips
		ST301	3,840 tips
Tubes	ScreenTape 0.2ml sample tubes and lids (in strips of 8 tubes)	ST221	250 strips
96 well plates	ScreenTape 96 well sample plate	ST219	25 plates

## Additional equipment recommended for optimal use of the Lab901 TapeStation

The table below details equipment that we recommend you use with the Lab901 TapeStation to achieve optimal performance and ease of use. All products listed below are available to purchase from Lab901 Ltd.

Equipment	Description	Catalogue Number
ScreenTape Micro-centrifuge (for 0.2ml tube strips, 220V)	For use during sample preparation to ensure that samples are free of bubbles and placed at the bottom of the sample tube	ST302
ScreenTape Micro-centrifuge (for 0.2ml tube strips, 110V)		ST305
ScreenTape Heating block (for 0.2 ml tube strips, 220V)	For use during sample preparation to allow optimal sample heating during staining and denaturing steps	ST303
ScreenTape Heating block (for 0.2 ml tube strips, 110V)		ST306
ScreenTape multichannel pipette	For loading ScreenTape filtered loading tips into the TapeStation	ST209
ScreenTape sample barcode scanner	To scan barcode labelled samples for efficient and reliable sample description entry	ST222