

ScreenTape P200 is a fully automated electrophoresis platform for proteins between 10 and 200kDa. This novel system allows you to analyse 16 protein samples in 12 minutes with no need for gel or buffer handling. The ScreenTape platform delivers accurate results for molecular weight sizing, protein purity and protein quantification. Complex protein purification workflows can be streamlined and accelerated by switching to Lab901's automated electrophoresis system for all your protein QC steps. With P200 the protein analysis bottleneck is eliminated and results are archived and can be exported to Word, Excel and a GLP standard report in a single mouse click.

